

Lifting Trolleys

Lifting (Scissor Lift) Trolleys

We manufacture a custom range of high quality Lifting Trolleys, that are built to exacting standards. Fully compatible with the ParMED One-Tray system, our trolleys can accommodate upto 750mm wide Body Trays. Fabricated in Stainless Steel, our custom design ensures operator ergonomics, smooth operation and a well engineered reliability.

Standard features include:

- Stainless Steel Construction
- Fully electric operation (hydraulic actuation is available as an option)
- Easy maintenance 24V Battery operation with local components

- Height Range of 400mm to 1780mm as standard

- Minimum Safe Working Load of 200kg
- Non-corrosive rollers in Galvanised Steel or HDPE
- 4 off Swivel Castors. 3 off with Brake and 1 off with Directional Locking

Optional features include:

- Powered Tray roller for Automatic tray retraction from racking
- Transaxle Wheel drive for motorised forward and reverse functions
- Upgraded Load rating of upto 300 kg SWL
- Customised Height Range
- Custom length profile for tight corners and rooms
- Load Cell weighing device for on trolley weight display